
**stypro®**  
Haemostypticum

Aufbewahrung trocken, zwischen +10°C und +25°C!  
Store dry, between +10°C and +25°C!  
Stocker dans un endroit sec entre +10°C et +25°C!  
iConservar en lugar seco, entre +10°C y +25°C!  
Conservare in un luogo asciutto a una temperatura!  
compresa tra +10°C e +25°C!  
Armazenar em local seco, entre +10°C and +25°C!

LOT

Chargenbezeichnung  
Batch number  
Mention de la charge  
Designación del lote  
Denominazione del lotto  
Número do lote

Verwendbar bis  
Expiry date  
A utiliser avant le  
Utilizable hasta  
Da utilizzare entro  
Data de Vencimento

Nur zum einmaligen Gebrauch!  
For single use only!  
Exclusivement pour usage unique!  
¡Solo para uso único!  
Solo monouso!  
Uso único

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Instructions for Use

GB

**stypro®** Sterile and haemostatic gelatine sponge for single use  
**stypro®** is sterilized by gamma radiations

Directions:

stypro® is applied directly in situ or with a slight pressure. stypro® being dry or moisturized previously with a sterile saline solution is applied pressing slightly on the clean wound, until the adhesion exerted by fibrin is obtained. The sponge can be left in situ, or removed just after the haemostasis has occurred. stypro® is completely absorbed within 3 weeks approximately. The amount of stypro® to be applied depends proportionally on the extent of the bleeding area. The sponge can be cut into strips of a desired size with sterile scissors, without occurrences of crumbings or breaks. stypro® has a poor adhesion to gauze pads or retractors.

Warnings:

Before use, verify the integrity of the package. Sterility is guaranteed if the package is maintained dry and not damaged. Apply stypro® directly after the opening after the package. Use stypro® before the expiration date printed on the package. When applying the product, strict aseptic conditions should be maintained. Strips of stypro® cannot be re-sterilized, and should be discarded. The haemostatic activity of stypro® however does not allow neglecting normal haemostatic procedures (suture, pressure and ligature). When stypro® is introduced into cavities it is recommended that the product is compressed, since the liquid absorption causes the haemostatic sponge to expand and this expansion may exert a pressure on the surrounding structures. In Neurosurgery the haemostatic sponge should not be left in closed cavities, since it may cause compression of the nervous tissue when it expands. stypro®, like other haemostatic gelatine sponges is intended to foster the stoppage of bleeding; any other use of the product must carefully be assessed in advance by the user. Keep out of reach of children.

Contraindications:

Not to be used on infected wounds

Storage:

Store in dry conditions

Preparation:

- stypro® Standard:** 80 x 50 x 10 mm: packages with 20 or 5 units
- stypro® Special:** 80 x 50 x 1 mm: packages with 20 or 5 units
- stypro® Tampon:** 80 x Ø 30 mm: packages with 20 or 2 units
- stypro® Cubus:** 10 x 10 x 10 mm: packages with 30 units
- stypro® Sheet:** 80 x 50 x 3 mm: packages with 20 units
- stypro® Strip:** 50 x 10 x 10 mm: packages with 5 units

The manufacturing company declines all responsibility for complaints deriving from an improper use of the product.

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